



ALS Prague Thermal Desorption

The purpose is to give an idea on how the thermal desorption process works, to pinpoint the key issues involving sample handling and the sampling process and to show method validation progress in the ALS Prague Laboratory



RIGHT SOLUTIONS | RIGHT PARTNER

List of topics




- Principle of thermal desorption
- Initial key steps of TO-17 method
- Sorbent selection process
- Importance of laboratory and field blanks
- Tube preparation and tube handling
- Summary of ALS method validation progress
- Sampling volume

Principle of 2 stage thermal desorption process



Method TO-17 5.2 Two-stage Thermal Desorption - the process of thermally desorbing analytes from a solid or liquid matrix, reconcentrating them on a focusing tube and then rapidly heating the tube to „inject” the concentrated compounds into the GC system in a narrow band of vapor compatible with high resolution capillary gas chromatography.





TO-17 2.3 Initial key steps of method

2.3.1 Selection of a sorbent or sorbent mix tailored for a target compound list, data quality objectives and sampling environment.

2.3.2 Screening the sampling location for VOCs by taking single tube samples to allow estimates of the nature and amount of sample gases.

2.3.3 Initial sampling sequences with two tubes at nominally 1 and 4 liter total sample volumes (or appropriate proportional scaling of these volumes to fit the target list and monitoring objectives).

[Note: EPA requires the use of distributed volume pairs for monitoring to insure high quality data. However, in situations where acceptable data have been routinely obtained through use of distributed volume pairs and the ambient air is considered well characterized, cost considerations may warrant single tube sampling. Any attendant risk to data quality objectives is the responsibility of the project's decision maker.]

- A successful analysis depends on the correct sorbent selection for the desired compound range.
- One of the illustrative criteria for sorbent selection is the challenging volume (the volume pulled through the tube needed to release an analyte)
- The following 3 slides demonstrate proposed composition of sorbents contained in a sampling tube which is able to capture all required analytes. Combination of multiple sorbents has several advantages over single sorbent tubes.

Challenge Volume (Liters)

	0.2	1	5	10	20	100
Halocarbon 12	Red	Red	Red	Red	Red	Red
Chloromethane	Red	Red	Red	Red	Red	Red
Halocarbon 114	Red	Red	Red	Red	Red	Red
Vinyl Chloride	Red	Red	Red	Red	Red	Red
1,3-Butadiene	Red	Red	Red	Red	Red	Red
Bromomethane	Red	Red	Red	Red	Red	Red
Ethyl Chloride	Red	Red	Red	Red	Red	Red
Halocarbon 11	Yellow	Red	Red	Red	Red	Red
Acrylonitrile	Green	Green	Yellow	Yellow	Yellow	Yellow
1,1-Dichloroethylene	Red	Red	Red	Red	Red	Red
Methylene Chloride	Green	Red	Red	Red	Red	Red
3-Chloropropylene	Green	Yellow	Red	Red	Red	Red
Halocarbon 113	Red	Red	Red	Red	Red	Red
1,1 Dichloroethane	Green	Green	Red	Red	Red	Red
cis-1,2 Dichloroethane	Green	Red	Red	Red	Red	Red
Chloroform	Green	Green	Red	Red	Red	Red
1,2 Dichloroethane	Green	Green	Yellow	Red	Red	Red
1,1,1 Trichloroethane	Green	Yellow	Yellow	Red	Red	Red
Benzene	Green	Green	Red	Red	Red	Red
Carbon Tetrachloride	Green	Yellow	Red	Red	Red	Red
1,2-Dichloropropane	Green	Green	Yellow	Yellow	Yellow	Red
Trichloroethylene	Green	Green	Yellow	Yellow	Red	Red
cis-1,3 Dichloropropene	Green	Green	Green	Green	Yellow	Red
trans-1,3-Dichloropropene	Green	Green	Green	Green	Green	Red
1,1,2-Trichloroethane	Green	Green	Green	Green	Yellow	Red
Toluene	Green	Green	Green	Green	Green	Red
1,2-Dibromoethane	Green	Green	Green	Green	Green	Red
Tetrachloroethylene	Green	Green	Green	Green	Yellow	Red
Chlorobenzene	Green	Green	Green	Green	Green	Yellow
Ethylbenzene	Green	Green	Green	Green	Green	Green
m,p-Xylene	Green	Green	Green	Green	Green	Green
Styrene	Green	Green	Green	Green	Green	Green
1,1,2,2-Tetrachlorethylene	Green	Green	Green	Green	Green	Green
o-xylene	Green	Green	Green	Green	Green	Green
4-Ethyltoluene	Green	Green	Green	Green	Green	Green
1,3,5-Trimethylbenzene	Green	Green	Green	Green	Green	Green
1,2,4-Trimethylbenzene	Green	Green	Green	Green	Green	Green
1,3-Dichlorobenzene	Green	Green	Green	Green	Green	Green
1,4-Dichlorobenzene	Green	Green	Green	Green	Green	Green
1,2-Dichlorobenzene	Green	Green	Green	Green	Green	Green
1,2,4-Trichlorobenzene	Green	Green	Green	Green	Green	Green
Hexachloro-1,3-butadiene	Green	Green	Green	Green	Green	Green

Tenax TA challenging volumes

Comment: Tenax TA is not suitable for all compounds from Chloroform and above

Performance Key

Safe to use: Recovery is greater than 80%
Caution: Recovery is between 21 to 79%
Not Recommended: Recovery is less than 20%

* indicates this analyte was strongly adsorbed



Carbograph1TD
challenging
volumes

	0.2	1	5	10	20	100
Halocarbon 12	Red	Red	Red	Red	Red	Red
Chloromethane	Red	Red	Red	Red	Red	Red
Halocarbon 114	Green	Red	Red	Red	Red	Red
Vinyl chloride	Red	Red	Red	Red	Red	Red
1,3-Butadiene	Green	Red	Red	Red	Red	Red
Bromomethane	Red	Red	Red	Red	Red	Red
Chloroethane	Red	Red	Red	Red	Red	Red
Halocarbon 11	Green	Red	Red	Red	Red	Red
Acrylonitrile	Green	Green	Green	Green	Red	Red
1,1-Dichloroethene	Green	Yellow	Red	Red	Red	Red
Methylene chloride	Red	Red	Red	Red	Red	Red
3-Chloropropene	Green	Yellow	Red	Red	Red	Red
Halocarbon 113	Green	Green	Green	Yellow	Red	Red
1,1-Dichloroethane	Green	Green	Red	Red	Red	Red
cis-1,2-Dichloroethene	Green	Green	Red	Red	Red	Red
Chloroform	Green	Green	Red	Red	Red	Red
1,2-Dichloroethane	Green	Green	Green	Red	Red	Red
1,1,1-Trichloroethane	Green	Green	Green	Green	Green	Red
Benzene	Green	Green	Green	Green	Green	Green
Carbon tetrachloride	Green	Green	Green	Green	Green	Red
1,2-Dichloropropane	Green	Green	Green	Green	Green	Red
Trichloroethene	Green	Green	Green	Green	Green	Yellow
cis-1,3-Dichloropropene	Green	Green	Green	Green	Green	Red
trans-1,3-Dichloropropene	Green	Green	Green	Green	Green	Red
1,1,2-Trichloroethane	Green	Green	Green	Green	Green	Green
Toluene	Green	Green	Green	Green	Green	Green
1,2-Dibromoethane	Green	Green	Green	Green	Green	Yellow
Tetrachloroethene	Green	Green	Green	Green	Green	Green
Chlorobenzene	Green	Green	Green	Green	Green	Green
Ethylbenzene	Green	Green	Green	Green	Green	Green
m & p-Xylene	Green	Green	Green	Green	Green	Green
Styrene	Green	Green	Green	Green	Green	Green
1,1,2,2-Tetrachlorethane	Green	Green	Green	Green	Green	Green
o-Xylene	Green	Green	Green	Green	Green	Green
4-Ethyltoluene	Green	Green	Green	Green	Green	Green
1,3,5-Trimethylbenzene	Green	Green	Green	Green	Green	Green
1,2,4-Trimethylbenzene	Green	Green	Green	Green	Green	Green
1,3-Dichlorobenzene	Green	Green	Green	Green	Green	Green
1,4-Dichlorobenzene	Green	Green	Green	Green	Green	Green
1,2-Dichlorobenzene	Green	Green	Green	Green	Green	Green
1,2,4-Trichlorobenzene	Green	Green	Green	Green	Green	Green
Hexachlorobutadiene	Green	Green	Green	Green	Green	Green

Performance Key

Safe to use: Recovery is greater than 80%
Caution: Recovery is between 21 to 79%
Not Recommended: Recovery is less than 20%

* indicates this analyte was strongly adsorbed



Carboxen 1003 challenging volumes

Challenge Volume (Liters)

	0.2	1	5	10	20	100
Halocarbon 12	Green	Green	Green	Green	Green	Red
Chloromethane	Green	Yellow	Red	Red	Red	Red
Halocarbon 114	Green	Green	Green	Green	Green	Green
Vinyl chloride	Green	Green	Green	Green	Green	Yellow
1,3-Butadiene	Yellow	Yellow	Yellow	Yellow	Yellow	Yellow
Bromomethane	Green	Yellow	Red	Red	Red	Red
Chloroethane	Green	Green	Green	Green	Green	Red
Halocarbon 11	Green	Green	Green	Green	Green	Green
Acrylonitrile	Green	Green	Green	Green	Green	Green
1,1-Dichloroethene	Green	Green	Green	Green	Green	Green
Methylene chloride	Green	Green	Green	Green	Green	Yellow
3-Chloropropene	Yellow	Yellow	Yellow	Yellow	Yellow	Yellow
Halocarbon 113	Green	Green	Green	Green	Green	Green
1,1-Dichloroethane	Green	Green	Green	Green	Green	Green
cis-1,2-Dichloroethene	Green	Green	Green	Green	Green	Green
Chloroform	Green	Green	Green	Green	Green	Yellow
1,2-Dichloroethane	Green	Green	Green	Green	Green	Green
1,1,1-Trichloroethane	Yellow	Yellow	Yellow	Yellow	Yellow	Yellow
Benzene *	Green	Green	Green	Green	Green	Yellow
Carbon tetrachloride	Yellow	Yellow	Yellow	Yellow	Yellow	Red
1,2-Dichloropropane	Yellow	Yellow	Yellow	Yellow	Yellow	Yellow
Trichloroethene	Green	Green	Green	Green	Green	Green
cis-1,3-Dichloropropene	Red	Red	Red	Red	Red	Red
trans-1,3-Dichloropropene	Red	Red	Red	Red	Red	Red
1,1,2-Trichloroethane	Yellow	Yellow	Yellow	Yellow	Yellow	Yellow
Toluene *	Red	Red	Red	Red	Red	Red
1,2-Dibromoethane	Red	Red	Red	Red	Red	Red
Tetrachloroethene *	Yellow	Yellow	Yellow	Yellow	Yellow	Yellow
Chlorobenzene *	Yellow	Yellow	Yellow	Yellow	Yellow	Yellow
Ethylbenzene *	Yellow	Yellow	Yellow	Yellow	Yellow	Red
m & p-Xylene *	Yellow	Yellow	Yellow	Yellow	Yellow	Red
Styrene *	Red	Red	Red	Red	Red	Red
1,1,2,2-Tetrachlorethane	Red	Red	Red	Red	Red	Red
o-Xylene *	Yellow	Yellow	Yellow	Yellow	Yellow	Yellow
4-Ethyltoluene	Red	Red	Red	Red	Red	Red
1,3,5-Trimethylbenzene	Red	Red	Red	Red	Red	Red
1,2,4-Trimethylbenzene	Red	Red	Red	Red	Red	Red
1,3-Dichlorobenzene *	Red	Red	Red	Red	Red	Red
1,4-Dichlorobenzene *	Red	Red	Red	Red	Red	Red
1,2-Dichlorobenzene *	Red	Red	Red	Red	Red	Red
1,2,4-Trichlorobenzene	Red	Red	Red	Red	Red	Red
Hexachlorobutadiene	Red	Red	Red	Red	Red	Red



Performance Key

Safe to use: Recovery is greater than 80%
Caution: Recovery is between 21 to 79%
Not Recommended: Recovery is less than 20%

* indicates this analyte was strongly adsorbed





**ALS Prague:
Comparison
between Tenax TA
and multisorbent
tubes – challenge
volumes tested:
0.2, 0.5, 1 and 2
Litres**

Compound name	TENAX				Multisorbent			
	0.2L	0.5L	1L	2L	0.2L	0.5L	1L	2L
1,2-Dichloro-1,1,2,2-tetrafluoroethane (CFC-114)	Red	Red	Red	Red	Green	Green	Green	Green
Chloromethane	Red	Red	Red	Red	Green	Green	Green	Green
Vinyl Chloride	Red	Red	Red	Red	Green	Green	Green	Green
1,3-Butadiene	Red	Red	Red	Red	Green	Green	Green	Green
Bromomethane	Red	Red	Red	Red	Green	Green	Green	Green
Chloroethane	Red	Red	Red	Red	Green	Green	Green	Green
Trichlorofluoromethane	Yellow	Red	Red	Red	Green	Green	Green	Green
Ethanol	Yellow	Red	Red	Red	Green	Green	Green	Green
1,1-Dichloroethene	Red	Red	Red	Red	Green	Green	Green	Green
1,1,2-Trichloro-1,2,2-trifluoroethane (CFC-113)	Red	Red	Red	Red	Green	Green	Green	Green
Acetone	Red	Red	Red	Red	Green	Green	Green	Green
Carbon Disulfide	Green	Green	Yellow	Yellow	Green	Green	Green	Green
2-Propanol (Isopropyl Alcohol)	Yellow	Red	Red	Red	Green	Green	Green	Green
Dichloromethane (Methylene Chloride)	Green	Green	Red	Red	Green	Green	Green	Green
trans-1,2-Dichloroethene	Green	Green	Yellow	Yellow	Green	Green	Green	Green
n-Hexane	Green	Green	Green	Green	Green	Green	Green	Green
1,1-Dichloroethane	Green	Green	Yellow	Yellow	Green	Green	Green	Green
2-Butanone (MEK)	Yellow	Red	Red	Red	Green	Green	Green	Green
cis-1,2-Dichloroethene	Green	Green	Yellow	Yellow	Green	Green	Green	Green
Tetrahydrofuran (THF)	Green	Green	Green	Green	Green	Green	Green	Green
Chloroform	Green	Green	Green	Green	Green	Green	Green	Green
1,1,1-Trichloroethane	Yellow	Red	Red	Red	Green	Green	Green	Green
Cyclohexane	Yellow	Red	Red	Red	Green	Green	Green	Green
Tetrachloromethane (Carbon Tetrachloride)	Yellow	Red	Red	Red	Green	Green	Green	Green
Benzene	Green	Green	Green	Green	Green	Green	Green	Green
1,2-Dichloroethane	Green	Green	Green	Green	Green	Green	Green	Green
2,2,4-Trimethylpentane (Isooctane)	Yellow	Red	Red	Red	Green	Green	Green	Green
n-Heptane	Green	Green	Green	Green	Green	Green	Green	Green
Trichloroethene	Green	Green	Green	Green	Green	Green	Green	Green
Methylcyclohexane	Yellow	Red	Red	Red	Green	Green	Green	Green
1,2-Dichloropropane	Green	Green	Green	Green	Green	Green	Green	Green
1,4-Dioxane	Green	Green	Green	Green	Green	Green	Green	Green
Methyl iso-Butyl Ketone (4-Methyl-2-pentanone)	Green	Green	Green	Green	Green	Green	Green	Green
Toluene	Green	Green	Green	Green	Green	Green	Green	Green
1,1,2-Trichloroethane	Green	Green	Green	Green	Green	Green	Green	Green
Tetrachloroethene	Green	Green	Green	Green	Green	Green	Green	Green
2-Hexanone	Yellow	Red	Red	Red	Green	Green	Green	Green
Chlorobenzene	Green	Green	Green	Green	Green	Green	Green	Green
Ethylbenzene	Green	Green	Green	Green	Green	Green	Green	Green
meta-Xylene & para-Xylene	Green	Green	Green	Green	Green	Green	Green	Green
o-Xylene	Green	Green	Green	Green	Green	Green	Green	Green
Styrene	Green	Green	Green	Green	Green	Green	Green	Green
Isopropylbenzene (Cumene)	Green	Green	Green	Green	Green	Green	Green	Green
1,1,2,2-Tetrachloroethane	Green	Green	Green	Green	Green	Green	Green	Green
n-Propylbenzene	Green	Green	Green	Green	Green	Green	Green	Green
4-Ethyltoluene	Green	Green	Green	Green	Green	Green	Green	Green
1,3,5-Trimethylbenzene	Green	Green	Green	Green	Green	Green	Green	Green
1,2,4-Trimethylbenzene	Green	Green	Green	Green	Green	Green	Green	Green
1,3-Dichlorobenzene	Green	Green	Green	Green	Green	Green	Green	Green
1,4-Dichlorobenzene	Green	Green	Green	Green	Green	Green	Green	Green
1,2-Dichlorobenzene	Green	Green	Green	Green	Green	Green	Green	Green
1,3,5-Trichlorobenzene	Green	Green	Green	Green	Green	Green	Green	Green
1,2,4-Trichlorobenzene	Green	Green	Green	Green	Green	Green	Green	Green
Hexachlorobutadiene	Green	Green	Green	Green	Green	Green	Green	Green
1,2,3-Trichlorobenzene	Green	Green	Green	Green	Green	Green	Green	Green

Less than 20% recovery
Between 20% and 80% recovery
Better than 80% recovery

Comment:
The first 5
analytes had
the recovery
less than 1%
on Tenax TA

Laboratory blanks (TO-17/Markes)



- Laboratory blanks are prepared in the same way and at the same time as tubes to be used for a monitoring exercise, but are never removed from the laboratory. They are simply kept inside a suitable storage container throughout the sample collection process and analysed with the sampled tubes. Laboratory blanks are used to check tube conditioning procedures and inherent sorbent stability.



Field blanks (TO-17/Markes)

- Field blanks are prepared as laboratory blanks above, but are transported with the tubes to be used for sample collection. Once at the monitoring location, all the processes carried out on the sampling tubes (e.g. unwrapping, uncapping, recapping, rewrapping, etc.), except sample collection itself, should also be carried out on the field blanks. They provide important controls for all tube handling, storage and transport procedures.

Tube preparation, handling and storage



- Tube preparation begins with an adequate conditioning process
- 10% of conditioned tubes within a conditioned batch are tested
- Tubes are capped with brass/PTFE storage caps using Cap Lok tools
- Tubes are stored in an air tight container



Tube preparation, handling and storage



- If the multi-sorbent tubes are to be stored for longer than one week, or will be exposed to low temperatures (an aircraft hold, by rail / road overnight during cold weather) then the following procedure must be undertaken:
- take the samples in the normal way and cap with ¼-inch brass storage caps in the field as usual
- place the tubes in an air tight container and place in a clean refrigerator
- when the tubes have reached the 'cold' temperature, briefly remove them from the refrigerator and check and re-tighten the caps (due to the difference in the thermal properties of brass and stainless steel the caps may have become slightly loose).
- return the samples, in their container, to the clean refrigerator

Calibration Range ALS Prague



- Calibration is done on multi-sorbent tube (TENAX TA, Carbograph 1TD, Carboxen 1003)
- Calibration range is 2 ng/tube - 1000 ng/tube
 - exception: Bromomethane, CFC-114 (upper limit 100ng/tube)
 - Exception: Chloromethan (upper limit 400ng/tube)

ALS Prague method limits of quantification – fortified blanks at 2ng/tube



Compound name	LOQ [ng/tube]	Compound name	LOQ [ng/tube]
1.2-Dichloro-1.1.2.2-tetrafluoroethane (CFC-114)	4.0	Trichloroethene	4.0
Chloromethane	14.0	Methylcyclohexane	4.0
Vinyl Chloride	2.0	1.2-Dichloropropane	4.0
1.3-Butadiene	2.0	1.4-Dioxane	10.0
Bromomethane	4.0	Methyl iso-Butyl Ketone	4.0
Chloroethane	4.0	Toluene	4.0
Trichlorofluoromethane	4.0	1.1.2-Trichloroethane	4.0
Ethanol	10.0	Tetrachloroethene	4.0
1.1-Dichloroethene	4.0	2-Hexanone	4.0
1.1.2-Trichloro-1.2.2-trifluoroethane (CFC-113)	4.0	Chlorobenzene	4.0
Acetone	20.0	Ethylbenzene	4.0
Carbon Disulfide	20.0	meta-Xylene a para-Xylene	4.0
2-Propanol (Isopropyl Alcohol)	4.0	o-Xylene	4.0
Dichloromethane (Methylene Chloride)	10.0	Styrene	4.0
trans-1,2-Dichloroethene	4.0	Isopropylbenzene (Cumene)	4.0
n-Hexane	4.0	1.1.2.2-Tetrachloroethane	4.0
1.1-Dichloroethane	4.0	n-Propylbenzene	4.0
2-Butanone (MEK)	20.0	4-Ethyltoluene	4.0
cis-1,2-Dichloroethene	4.0	1.3.5-Trimethylbenzene	4.0
Tetrahydrofuran (THF)	4.0	1.2.4-Trimethylbenzene	4.0
Chloroform	4.0	1.3-Dichlorobenzene	4.0
1.1.1-Trichloroethane	4.0	1.4-Dichlorobenzene	4.0
Cyclohexane	2.0	1.2-Dichlorobenzene	4.0
Tetrachloromethane (Carbon Tetrachloride)	4.0	1.3.5-Trichlorobenzene	10.0
Benzene	4.0	1.2.4-Trichlorobenzene	10.0
1.2-Dichloroethane	4.0	Hexachlorobutadiene	10.0
2.2.4-Trimethylpentane (Isooctane)	4.0	1.2.3-Trichlorobenzene	10.0
n-Heptane	4.0		

ALS Prague repeatability – fortified blanks at 2ng/tube



Compound name	RSD	BIAS	Compound name	RSD	BIAS
1.2-Dichloro-1.1.2.2-tetrafluoroethane (CFC-114)	8.4%	10.9%	Trichloroethene	4.0%	2.4%
Chloromethane	31.5%	23.9%	Methylcyclohexane	3.3%	8.4%
Vinyl Chloride	4.3%	0.2%	1.2-Dichloropropane	4.7%	2.8%
1.3-Butadiene	11.6%	27.2%	1.4-Dioxane	6.6%	6.3%
Bromomethane	4.1%	0.3%	Methyl iso-Butyl Ketone (4-Methyl-2-pentanone)	4.6%	1.5%
Chloroethane	5.2%	4.1%	Toluene	5.3%	6.1%
Trichlorofluoromethane	6.4%	4.9%	1.1.2-Trichloroethane	4.8%	14.7%
Ethanol	31.7%	31.5%	Tetrachloroethene	13.7%	20.6%
1.1-Dichloroethene	6.9%	4.2%	2-Hexanone	3.8%	7.6%
1.1.2-Trichloro-1.2.2-trifluoroethane (CFC-113)	5.3%	1.1%	Chlorobenzene	2.7%	12.2%
Acetone	30.2%	42.4%	Ethylbenzene	4.0%	11.1%
Carbon Disulfide	3.0%	6.7%	meta-Xylene a para-Xylene	4.9%	11.5%
2-Propanol (Isopropyl Alcohol)	10.5%	3.6%	o-Xylene	4.7%	24.7%
Dichloromethane (Methylene Chloride)	6.5%	4.2%	Styrene	3.5%	14.9%
trans-1,2-Dichloroethene	5.2%	1.3%	Isopropylbenzene (Cumene)	6.1%	0.8%
n-Hexane	5.5%	7.5%	1.1.2.2-Tetrachloroethane	4.7%	16.7%
1.1-Dichloroethane	6.9%	4.0%	n-Propylbenzene	13.8%	38.7%
2-Butanone (MEK)	5.0%	3.4%	4-Ethyltoluene	4.8%	17.7%
cis-1,2-Dichloroethene	12.6%	14.0%	1.3.5-Trimethylbenzene	7.7%	15.7%
Tetrahydrofuran (THF)	5.3%	1.6%	1.2.4-Trimethylbenzene	4.5%	11.4%
Chloroform	5.7%	7.3%	1.3-Dichlorobenzene	6.5%	22.0%
1.1.1-Trichloroethane	5.6%	5.8%	1.4-Dichlorobenzene	4.2%	4.0%
Cyclohexane	4.8%	1.1%	1.2-Dichlorobenzene	5.8%	28.2%
Tetrachloromethane (Carbon Tetrachloride)	5.4%	3.8%	1.3.5-Trichlorobenzene	5.3%	38.9%
Benzene	7.0%	4.3%	1.2.4-Trichlorobenzene	5.5%	32.9%
1.2-Dichloroethane	7.5%	4.6%	Hexachlorobutadiene	4.1%	38.6%
2.2.4-Trimethylpentane (Isooctane)	4.8%	2.0%	1.2.3-Trichlorobenzene	6.5%	32.5%
n-Heptane	5.9%	7.3%			

ALS Prague repeatability – fortified blanks at 100ng/tube



Compound name	RSD	BIAS	Compound name	RSD	BIAS
1.2-Dichloro-1.1.2.2-tetrafluoroethane (CFC-114)	6.1%	4.3%	Trichloroethene	2.0%	7.2%
Chloromethane	7.4%	14.7%	Methylcyclohexane	2.3%	3.4%
Vinyl Chloride	5.7%	11.9%	1.2-Dichloropropane	2.7%	0.7%
1.3-Butadiene	9.0%	5.8%	1.4-Dioxane	2.3%	6.4%
Bromomethane	7.8%	1.2%	Methyl iso-Butyl Ketone (4-Methyl-2-pentanone)	2.4%	3.2%
Chloroethane	5.9%	8.2%	Toluene	2.0%	9.9%
Trichlorofluoromethane	4.2%	7.5%	1.1.2-Trichloroethane	1.9%	5.0%
Ethanol	3.7%	4.1%	Tetrachloroethene	1.8%	8.5%
1.1-Dichloroethene	4.0%	6.8%	2-Hexanone	5.3%	10.1%
1.1.2-Trichloro-1.2.2-trifluoroethane (CFC-113)	2.8%	1.9%	Chlorobenzene	1.6%	7.1%
Acetone	4.6%	11.0%	Ethylbenzene	2.1%	5.0%
Carbon Disulfide	3.4%	7.3%	meta-Xylene a para-Xylene	2.4%	3.8%
2-Propanol (Isopropyl Alcohol)	4.7%	6.1%	o-Xylene	1.9%	6.6%
Dichloromethane (Methylene Chloride)	3.2%	0.4%	Styrene	2.1%	9.3%
trans-1,2-Dichloroethene	2.6%	2.3%	Isopropylbenzene (Cumene)	2.4%	7.8%
n-Hexane	2.9%	5.7%	1.1.2.2-Tetrachloroethane	2.1%	5.2%
1.1-Dichloroethane	2.7%	1.1%	n-Propylbenzene	2.4%	6.4%
2-Butanone (MEK)	3.6%	6.9%	4-Ethyltoluene	2.4%	10.2%
cis-1,2-Dichloroethene	2.0%	3.8%	1.3.5-Trimethylbenzene	2.3%	9.8%
Tetrahydrofuran (THF)	2.9%	4.9%	1.2.4-Trimethylbenzene	2.5%	6.9%
Chloroform	2.1%	1.4%	1.3-Dichlorobenzene	2.8%	9.2%
1.1.1-Trichloroethane	2.8%	0.6%	1.4-Dichlorobenzene	2.9%	10.6%
Cyclohexane	3.5%	4.0%	1.2-Dichlorobenzene	3.1%	9.0%
Tetrachloromethane (Carbon Tetrachloride)	2.7%	0.2%	1.3.5-Trichlorobenzene	2.3%	12.2%
Benzene	2.1%	2.6%	1.2.4-Trichlorobenzene	3.5%	16.2%
1.2-Dichloroethane	3.0%	4.0%	Hexachlorobutadiene	2.1%	12.1%
2.2.4-Trimethylpentane (Isooctane)	2.5%	0.3%	1.2.3-Trichlorobenzene	3.3%	13.2%
n-Heptane	3.1%	4.9%			

- What are the legislation limits for VOC' in Israel?
- How quickly needs to be sample taken?
- How many pumps are available for single sampling?
 - (Is duplicate determination possible?)
- Would it be possible to send out tubes to the sampling site in an air tight containers? Ideally 3 tubes per sampling site.
- Can we distribute Cap Lok tools with every sample?
- How long does it take to receive a sample in ALS Prague from the day the sample is taken?