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- Method is developed for determination non-polar extractable substances in soil, ground water, drinking water, waste water and leachate.
- Extraction solvent is tetrachloroethylene
- Co-extracting polar compounds are adsorbed from dried extract by shaking with silica gel
- IR spectrum of dried extract is measured in 10 mm cuvette with quartz sides in wavenumber range 2600–3300 cm^{-1} single beam IR instrument. The IR is connected with PC and the special software save the absorb spectrum of measured extraction solvent.

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- **Sample homogenization** – before weighing the sample is manually homogenized. In case of large particles presented in the sample (more than 5 mm) the sample is crushed in jaw crusher.
- **Sample weighing** – standard weigh of sample is 10 g. If it is obvious that the sample is contaminated, the sample weigh is lower.
- **Sample drying** – The sodium sulphate Na_2SO_4 is added to the sample and mixed. Dried sample should strew and do not stick on walls of extraction glass bottle.
- **Sample extraction** – 30 ml of tetrachloroethylene C_2Cl_4 is added to the dried sample and the mixture is shaken on horizontal shaker (speed 150 min^{-1}).
- **Extract separation** – extract is filtered thru paper filter. Filtrate is collected to the volumetric flask.
- **Extract clean-up** – is the same as water extract clean-up. Adding of activated silica gel and shaking.

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
- The concentration is calculated from maximum absorption values A_1 , A_2 and A_3 of the extract on three different wavenumbers, that correspond to vibration of C-H bound in characteristic atom groups:

A_1 (aromatics and olefinics =CH-)	$3\ 055 \pm 25\text{ cm}^{-1}$
A_2 (asymetrics -CH ₃)	$2\ 960 \pm 15\text{ cm}^{-1}$
A_3 (asymetrics -CH ₂)	$2\ 925 \pm 15\text{ cm}^{-1}$

- A lot of non-polar organic compounds can be determined by IR spectrometry. The sensitivity for aromatic compounds is smaller other techniques. There is no possibility to distinguish nature products by IR spectrometry. IR allows partially describe what type of hydrocarbons are presented in the sample. This information is based on absorption band ratios of -CH₃, -CH₂, =CH- groups. For better description of contamination is necessary to make GC analysis.

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- W-TPHW-IR

Wastewater

- required volume of sample 0,5 L
- LOQ 0.1 mg/L, measurement uncertainty (MU) 20 %

- W-TPH-IR

Drinking and groundwater

- required volume of sample 1 L
- LOQ 0.05 mg/L, MU 20 %

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Gasoline range (GRO)	>C6-C10
Diesel range (DRO)	>C11-C28
Lube oil range	>C28-C35

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
The GC/FID method is used for analysis of extractable hydrocarbons which are usually products of crude oil distillation. Refined fractions obtained during petroleum processing are mentioned in the following table.

	Boiling point (°C)	n-Alkanes
Gas hydrocarbons	< 5	C1 - C4
Light gasoline	30 - 85	C5 - C6
Heavy gasoline	85 - 180	C7 - C10
Paraffin oil	180 - 270	C11 - C15
Gas oil	270 - 370	C16 - C22
Vacuum spirits	370 - 550	C23 - C45
Vacuum waste	over 550	> C46

The GC/FID method is capable to analyze hydrocarbons in range of boiling point 85 to 350 °C. Low boiling fractions like Gas hydrocarbons or Light gasoline cannot be determined using by method GC/FID as well as heavy hydrocarbons like Vacuum Spirits and Vacuum Waste.

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
Fig. 1 Chromatogram of Standard Alkanes Mix 12 (Mixture of even Hydrocarbons from C10 to C40)

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